



Live and Dead Cell Double Staining Kit

Cat #: KTA1001

Size: 100 T/500 T/2000 T

| | | | |
|---|---|------------|--------------------------------------|
|  | Live and Dead Cell Double Staining Kit | | |
| REF | Cat #: KTA1001 | LOT | Lot #: Refer to product label |
| | Applicable samples: Cells | | |
|  | Storage: Stored at -20°C for 12 months, protected from light | | |

Assay Principle

Cell-mediated cytotoxicity is an important phenomenon characterized by cytolysis of a compromised cell in the body by immune system. Distinguishing between live and dead cells is very important for investigation of growth control and cell death. The Live/Dead Cell Double Staining Kit provides a convenient assay to evaluate the viability of cells, based on the simultaneous determination of live and dead cells with two probes that measure recognized parameters of cell health: plasma membrane integrity and intracellular esterase activity. The kit utilizes Ca-AM, a cell-permeable green fluorescent dye (Ex/Em=488/530 nm), to stain live cells and PI, a cell non-permeable red fluorescent dye (Ex/Em=535/617 nm), to stain dead cells.

Materials Supplied and Storage Conditions

| Kit components | Size | | | Storage conditions |
|--------------------|-------|--------|--------|-----------------------------|
| | 100 T | 500 T | 2000 T | |
| Ca-AM | 50 µL | 250 µL | 1 mL | -20°C, protected from light |
| PI | 50 µL | 250 µL | 1 mL | -20°C, protected from light |
| Assay Buffer (10×) | 5 mL | 25 mL | 100 mL | 4°C |

Materials Required but Not Supplied

- Microcentrifuge
- 24 well plate, Precision Pipettes, Disposable Pipette Tips
- Fluorescence Microscopy or Flow Cytometer
- Deionized Water, PBS

Reagent Preparation

Ca-AM: Keep on ice while using. Protect from light.

PI: Keep on ice while using. Protect from light.

1×Assay Buffer: Prepare 1×Assay Buffer by dilute 10× Assay Buffer with Deionized Water. Warm to 37°C before use.

Staining Solution: Mix 1 μ L Ca-AM and 1 μ L PI in each 1 mL Assay Buffer. Scale up accordingly for larger numbers of assays.

Assay Procedure

A. Quantification by Flow Cytometry

1. Treat cells with the desired method.

Note: We recommend keeping unstained control cells (i.e. without Ca-AM or PI staining) suspended in 1 \times Assay Buffer for both treated and untreated samples to set up the flow cytometer instrument.

2. For non-adherent cells, Collect 1-5 \times 10⁵ cells by centrifugation (4 $^{\circ}$ C, 300 g, 5 min). Wash with ice-cold PBS twice and discard the PBS. For adherent cells, using Trypsin (EDTA free) to digest cells firstly and then centrifugation.

3. Resuspend the cells pellet in 500 μ L Staining Solution.

4. Incubate the cells at 37 $^{\circ}$ C for 15-30 min in the dark.

5. Analyze the cells by flow cytometry.

B. Detection by Fluorescence Microscopy

1. For suspension cells: Follow the protocol for flow cytometry from step A.1 to step A.4 and place the cell suspension from Step A.4 on a glass slide. Cover the cells with a glass coverslip. Analyze cells by fluorescence microscopy using the appropriate filters as soon as possible.

2. For adherent cells: the suggested protocol is as below.

(1) Grow cells directly on a coverslip in 24 well plate. Incubate in a CO₂ Incubator at 37 $^{\circ}$ C for at least 24 h before treatment.

(2) Treat cells with the desired method. Prepare a negative control by incubating cells in the absence of inducing agent.

(3) Wash cells with PBS twice.

(4) Add 0.5 mL of Staining solution to cells and incubate at 37 $^{\circ}$ C for 15-30 min in the dark.

(5) Wash cells with PBS twice.

(6) Invert coverslip on a glass slide and visualize cells fluorescence microscopy using the appropriate filters as soon as possible.

Note: PI is a potential mutagen. Use appropriate precautions when handling this reagent.

Recommended Products

| Catalog No. | Product Name |
|-------------|---|
| KTA2020 | Cell Cycle Staining Kit |
| KTA1002 | Live Cell Tracking Kit (Green Fluorescence) |

Disclaimer

The reagent is only used in the field of scientific research, not suitable for clinical diagnosis or other purposes.