


Rayto RT-6900 ELISA Reader

電腦軟體快速操作指南



電腦系統前置設定



設定



搜尋裝置管理員 (Devices Manager)

設定

我的 Microsoft 帳戶

OneDrive 管理

Windows Update
上次檢查時間：31 分鐘前

裝置管理員

裝置管理員

系統
顯示器、音效、通知、電源

裝置
藍牙、印表機、滑鼠

電話
連結您的 Android、iPhone

網路和網際網路
Wi-Fi、飛航模式、VPN

個人化
背景、鎖定畫面、色彩

應用程式
解除安裝、預設值、選用功能

帳戶
您的帳戶、電子郵件、同步設定、工作、家庭

時間與語言
語音、地區、日期


遊戲
Xbox Game Bar、擷取、遊戲模式

輕鬆存取
朗讀程式、放大鏡、高對比

搜尋
尋找我的檔案、權限


隱私權
位置、相機、麥克風

1. 進入裝置管理員
2. 找到【連接埠 (COM和LPT)】並展開
3. 找到Prolific USB-to-Serial Comm Port (COM數字), 如**非COM 1**, 則按右鍵進入內容。




1. 選擇【連接埠設定】

2. 點入【進階】



3. 在【COM連接埠號碼】打開下拉選單，點選COM1。


4. 按【確定】完成設定




儀器前置設定

儀器開機後，請依下列指示操作，將儀器設定成**透過電腦操作的模式**：


1. 選擇【System Set】




2. 選擇【Sys. Set】



3. 選擇【PC Control】




4. 出現以下提示訊息即代表設定完成 (請勿按Cancel退出)



1. 啓動軟體

1.1 啟動軟體Rayto Microplate Reader Software



1.2 登入程式

Rayto Microplate Reader Software



UserName : Admin

Password : 888888

A screenshot of the software's login dialog box. The dialog has a title bar that says 'Login'. It contains two input fields: 'UserName' and 'Password'. The 'UserName' field is a dropdown menu with 'Admin' selected and highlighted in blue. The 'Password' field contains six black dots. At the bottom of the dialog are two buttons: 'Login' and 'Cancel'. A hand cursor is pointing at the 'Login' button.

Login

UserName Admin

Password ●●●●●●

Login Cancel

1.3 機器self-checking中



Rayto Microplate Reader Software

Filter wheel moves to wavelength 492.....



Item parameters

1.4 進入主畫面 (Test介面)Plate Direction Row ColumnStart sampleNo Well Position --

Plate current well position item information

Test parameters

Test template Testing process

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

2. 測吸光值

Item setting

QC setting

2.1 切換到Parameter介面，先設定Program。

Program

Full name

Unit

wavelength

Reference

Reagent

Setting

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Standard

Concentration

 Double standard

Qualitative way

Negative <

X

Positive >=

X

Setting

2.2 按New以新增Program

New

Save

Delete

Item setting

QC setting

2.3 輸入Program名稱

Program

Test 450

Full name

Unit

wavelength

405

None

Reference

Reagent

Setting

Measurement

End point meth

Calculation

ABS

Samples

1

Intervals time

Detection

Factor

Standard setting

Standard number

0

Double standard

Standard

Concentration

Standard	Concentration

Qualitative way

None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

Item setting

QC setting

2.4 選擇濃度單位 (僅測吸光值可不選)

Program Test 450

Full name

Unit

wavelength 405 None

Reference

Reagent U/L Setting

Measurement End point meth

Calculation ABS

Samples s/co

Intervals time

Detection

Factor ng/ml IU/ml

Standard setting

Standard number 0

Standard Concentration

Double standard

Qualitative way None

Negative < X

Positive > = X

Setting

New

Save

Delete

Item setting

QC setting

2.5 選擇波長一

Program Test 450

wavelength 450 None

Measurement 450

Intervals time 630

Full name

Reference

Calculation ABS

Detection

Unit

Reagent Setting

Samples 1

Factor

Standard setting

Standard number 0

Double standard

Standard Concentration

Qualitative way None

Negative < X

Positive >= X

Setting

New

Save

Delete

Item setting

QC setting

2.6 選擇波長二 (若無波長二則選None)

Program Test 450

wavelength 450

None

Measurement End point

Intervals time

None

405

450

492

630

Full name

Reference

Calculation ABS

Detection

Unit

Reagent

Samples 1

Factor

Standard setting

Standard number 0

Double standard

Standard

Concentration

Qualitative way

None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

Item setting

QC setting

2.7 輸入Reference：最小值0 -- 最大值1~4

Program Test 450

wavelength 450 None

Measurement End point meth

Intervals time

Full name

Reference

Calculation ABS

Detection

Unit

Reagent

Samples 1

Factor

Standard setting

Standard number 0

 Double standard

Standard

Concentration

Qualitative way None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

2.8 選擇吸光度測定方法 (一般選擇End point)

Item setting

QC setting

Program Test 450

wavelength 450 None

Measurement End point meth

Intervals time

Standard setting

Standard number 0

Double standard

Full name

Unit

Reference 0 -- 1

Reagent Setting

Calculation ABS

Samples 1

Detection

Factor

Qualitative way None

Negative < X

Positive > = X

Setting

New

Save

Delete

Standard	Concentration

Item setting

QC setting

2.9 選擇測定結果計算方法 (僅測吸光值選擇ABS)

Program Test 450

Full name

Unit

wavelength 450 None

Reference 0 -- 1

Reagent

Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Standard setting

Standard number 0

Standard

 Double standard

- ABS
- CutOff
- One point calibration
- %ABS
- Non-linear regression
- Linear regression
- Log-log regression
- Power regression
- Exponential regression
- Percent logarithm
- Cubic equation regression
- 4 parameters regression
- Quadratic regression

Qualitative way None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

Item setting

QC setting

2.10 選擇樣本重複數

Program Test 450

wavelength 450 None

Measurement End point meth

Intervals time

Full name

Reference 0 -- 1

Calculation ABS

Detection

Unit

Reagent Setting

Samples 3

Factor

- 1
- 2
- 3
- 4

Qualitative way None

Negative < X

Positive > = X

Setting

Standard setting

Standard number 0

Double standard

Standard	Concentration

New

Save

Delete

Item setting

QC setting

Program Test 450

Full name

Unit

wavelength 450 None

Reference 0 -- 1

Reagent Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Factor

Standard setting

Standard number 0

Double standard

Standard	Concentration

2.11 選擇定性方式

(僅測吸光非定性分析選擇None)

Qualitative way

- None
- None
- Positive threshold
- Reverse threshold

Negative <

Positive > = X

Setting

New

Save

Delete

Item setting

QC setting

Program Test 450

Full name

Unit

wavelength 450

None

Reference 0

1

Reagent

Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Factor

Standard setting

Standard number 0

 Double standard

Standard

Concentration

Qualitative way

None

Negative <

X

Positive >=

X

Setting

2.12 Program設定完成，按Save儲存。

New

Save

Delete



Item setting

QC setting

Test 450

Program Test 450

Full name

Unit

Reagent

Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Factor

Standard setting

Standard number 0

 Double standard

Standard

Concentration

Qualitative way

None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

2.13 儲存成功的Program會顯示在左側列表

Item parameters

2.14 切換到Test介面

Test 450

STD 630

**2.15 點選要套用的Program**Plate Direction Row ColumnStart sampleNo

Well Position

A

1

--

H

12

Plate current well position item information

Test parameters

Test template

Testing process

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Item parameters

Test 450

STD 630

2.16 點選Sample

2.17 選擇Sample排列方向

Plate Direction Row Column

Start sampleNo 0001

Well Position A 1 -- H 12 Setting

Plate current well position item information

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

2.18 點選空格排列Sample位置

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0010										
B	0002	0009										
C	0003											
D	0004											
E	0005											
F	0006											
G	0007											
H	0008											

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450 STD 630

2.19 快速設定Sample layout

- ① 選擇Sample排列方向
- ② 選擇Well起迄範圍
- ③ 按Setting
- ④ Sample layout套用完成

Plate Direction Row **①** Column

Start sampleNo 0001 **②**

Well Position A 1 -- H 2 **③** Setting

Plate current well position item informat

- 1
- 2**
- 3
- 4
- 5
- 6
- 7
- 8
- 9
- 10
- 11
- 12



Test parameters

Test template Setting Testing process default Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

④	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0009										
B	0002	0010										
C	0003	0011										
D	0004	0012										
E	0005	0013										
F	0006	0014										
G	0007	0015										
H	0008	0016										

Plate_Pop	Plate_Reset	Stop testing	Clear All	Result	ABS	Quantitative	Qualitative
Save	Print	Export	Send				

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo

Well Position --

Plate current well position item information

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

先按Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0010										
B	0002	0009										
C	0003											
D	0004											
E	0005											
F	0006											
G	0007											
H	0008											

2.20 如果Sample位置排錯了如何修正?

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

12

Setting

Plate current well position item information

Empty text area for plate current well position item information.

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001											
B	0002											
C	0003											
D	0004											
E	0005											
F	0006											
G	0007											
H	0008											



② 再點選排錯的格子即可完成清除

Item parameters

Test 450 STD 630

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

QC1

QC2

QC3

QC4

QC5

Blank

Clear

2.21 儲存Plate Template (可省略)

- ① 按Setting進入
- ② 按New
- ③ 自行命名Template
- ④ 按Save儲存

Plate Direction Row Column

Start sampleNo 0001

Well Position A 1 -- H 2 Setting

Plate current well position item information

Test template setting

Test template setting

Template 20200916-1

Template 20200916-1

New

Save

Delete

Back

New

Save

Delete

Back

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

2.22 日後可選擇已儲存的Plate template套用 (可省略)

Item parameters

Test 450

STD 630

Test template

20200916

Setting

Testing process

default

Setting

Sample

NC

20200916-1

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										
H	0015	0016										

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

2.23 設定shaking (可省略)

Item parameters

Test parameters

Test 450

STD 630

Test template 20200916

Setting

Testing process default

Setting



Testing process set

Testing process

default

Process settings

Delay 0 S

Shake Medium speed 5 S

Test

Step



Testing process set

Testing process

default

Process settings

Delay 0 S

Shake Medium speed 5 S

Test

Step

Shake 5s




Process name Shake 5s

New

Save

Delete

Back



Process name Shake 5s

New

Save

Delete

Back



Process name Shake 5s

New

Save

Delete

Back



2.24 可選擇shaking程序 (可省略)

Item parameters

Test 450

STD 630

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

default
Shake 5s

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										
H	0015	0016										

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo

Well Position

A

1

--

H

2

Setting

Plate current well position item information

2.25 按Plate_Reset 兩次

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										
H	0015	0016										

Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

2.26 按Plate_Pop, 托架彈出, 放入盤子 (A1朝右邊內側)。

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										
H	0015	0016										

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

2.27 按Test開始讀盤

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										
H	0015	0016										

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450

STD 630

Test parameters

Test template 20200916

Setting

Testing process Shake 5s

Setting

Sample NC PC PC2 CR QC1 QC2 QC3 QC4 QC5

STD Blank Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001	0002										
B	0003	0004										
C	0005	0006										
D	0007	0008										
E	0009	0010										
F	0011	0012										
G	0013	0014										

Plate Direction Row Column

Start sampleNo 0001

Well Position A 1 -- H 2 Setting

Plate current well position item information

(測讀過程中若要中止可按Stop testing)

Plate_Pop Plate_Reset Test **Stop testing** Clear All Result ABS Quantitative Qualitative
 Save Print Print Preview Export Send

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

Empty text area for plate current well position item information.

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0.000	0.000										
B	0.000	0.000										
C	0.000	0.000										
D	0.000	0.000										
E	0.000	0.001										
F	0.000	0.000										
G	0.000	0.000										
H	0.000	0.000										

2.28 測定完成，顯示結果。

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position A 1 -- H 2 Setting

Plate current well position item information

Test parameters

Test template 20200916 Setting Testing process Shake 5s Setting

Sample NC PC PC2 CR QC1 QC2 QC3 QC4 QC5

STD Blank Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001 0.000	0002 0.000										
B	0003 0.000	0004 0.000										
C	0005 0.000	0006 0.000										
D	0007 0.000	0008 0.000										
E	0009 0.000	0010 0.001										
F	0011 0.000	0012 0.000										
G	0013 0.000	0014 0.000										
H	0015 0.000	0016 0.000										

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

2.29 按 ABS顯示吸光值



Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

2.30 按Save進行存檔

(一定要先存檔才能輸出CSV檔)

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

1

2

3

4

5

6

7

8

9

10

11

12

A

0.000

0.000

B

C

D

E

F

0.000

0.000

G

0.000

0.000

H

0.000


0.000

Set Board Number

Plate 202009160001

Save

1



另存新檔

<< Program Files > MB > 200916

搜尋 200916

組合管理 新增資料夾

名稱 修改日期 類型 大小

沒有符合搜尋條件的項目

檔案名稱(N): 202009160001_200916_155030

存檔類型(T): csv file (*.csv)

存檔(S) 取消

20200916 Setting Testing process Shake 5s Setting

PC PC2 CR QC1 QC2 QC3 QC4 QC5

Clear

	3	4	5	6	7	8	9	10	11	12
F	0.000	0.000								
G	0.000	0.000								
H	0.000	0.000								

Well Position A 1 -- H 2 Setting

Plate current well position item information

Plate_Pop Plate_Reset Test Stop testing Clear All Result ABS Quantitative Qualitative

Save Print Print Preview Export Send

2.3 | 按Export輸出CSV檔

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

Empty text box for plate current well position item information.

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A	0.000	0.000										
B	0.000	0.000										
C	0.000	0.000										
D	0.000	0.000										
E	0.000	0.001										
F	0.000	0.000										
G	0.000	0.000										
H	0.000	0.000										



2.32 按Print Preview預覽列印PDF檔

2.33 按列印輸出成PDF檔

列印

Item name: Test 450 Full name: Method:ABS
 Range:0.000-1.000 Wavelength:450nm Reagent name:
 Manufacturer: Reagent Lot: Reagent Period:
 QC: QC Lot: QC Period:
 QCcontent: Temp/Humidity:°C %

Each square contains at most four rows, which represent the sample number, test results, quantitative results, and qualitative results.

	1	2	3	4	5	6	7	8	9	10	11	12
A	0001 0.000 0.000	0002 0.000 0.000										
B	0003 0.000 0.000	0004 0.000 0.000										
C	0005 0.000 0.000	0006 0.000 0.000										
D	0007 0.000 0.000	0008 0.000 0.000										
E	0009 0.000 0.000	0010 0.001 0.000										
F	0011 0.000 0.000	0012 0.000 0.000										
G	0013 0.000 0.000	0014 0.000 0.000										
H	0015 0.000 0.000	0016 0.000 0.000										

Print date:2020-09-16 Test time:2020-09-16 15:50:12 Tester:Admin Approver:

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo 0001

Well Position

A

1

--

H

2

Setting

Plate current well position item information

Test parameters

Test template

20200916

Setting

Testing process

Shake 5s

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

S

2.34 按Clear All將Plate內容全部清除



Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo

Well Position

--

Plate current well position item information

Test parameters

Test template

Testing process

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												



2.35 按Plate_Reset重設下一盤

3. Standard定量

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

Standard	Concentration

Qualitative way

Negative < X

Positive >= X

3.1 按New新增Program

New

Save

Delete



3.2 輸入Program名稱

Item setting

QC setting

Test 450

Program **STD 630**

Full name

Unit

wavelength 405 None

Reference

Reagent Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Factor

Standard setting

Standard number 0

Double standard

Standard	Concentration

Qualitative way None

Negative < X

Positive > = X

Setting

New

Save

Delete

3.3 選擇濃度單位 (定量必選)

Item setting

QC setting

Test 450

Program

Full name

wavelength

Reference --

Measurement

Calculation

Intervals time

Detection

Standard setting

Standard number

Double standard

Standard	Concentration

Unit

Reagent

Samples

Factor

Qualitative way

Negative < X

Positive > = X

Setting

New

Save

Delete

Item setting

QC setting

3.4 選擇波長一

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

Standard	Concentration

Qualitative way

Negative < X

Positive > = X

Item setting

QC setting

3.5 選擇波長二 (若無波長二則選None)

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

- None
- 405
- 450
- 492
- 630

Standard setting

Standard number

Double standard

Standard	Concentration

Qualitative way

Negative < X

Positive >= X

Item setting

QC setting

3.6 輸入Reference：最小值0 -- 最大值1~4

Test 450

Program Full name Unit wavelength Reference -- Reagent Measurement Calculation Samples Intervals time Detection Factor

Standard setting

Standard number Double standard

Standard	Concentration
----------	---------------

Qualitative way Negative < X Positive > = X

Item setting

QC setting

Test 450

Program STD 630

wavelength 630 405

Measurement End point meth

Intervals time

Standard setting

Standard number 0

 Double standard

New

3.7 如有需要紀錄使用的試劑資訊請按Setting (可省略)

Full name

Unit ng/ml

Reference 0 -- 1

Reagent

Calculation ABS

Samples 1

Setting

Reagent information setting

Reagent	Production lot	Manufacturer	Production Date	Reagent validity

Reagent

Production lot

Manufacturer

Production Date 2020/ 9/16

Reagent validity 2020/ 9/16

New

Save

Delete

Back

Item setting

QC setting

3.8 選擇吸光度測定方法（一般選擇End point）

Test 450

Program STD 630

Full name

Unit ng/ml

wavelength 630 405

Reference 0 -- 1

Reagent

Setting

Measurement End point meth

Calculation ABS

Samples 1

Intervals time

Detection

Factor

- End point method
- Two-point method
- Kinetic method

Standard setting

Standard number 0

 Double standard

Standard

Concentration

Qualitative way

None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

3.9 選擇測定結果計算方法

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Standard setting

Standard number

- Standard
- ABS
- CutOff
- One point calibration
- %ABS
- Non-linear regression
- Linear regression
- Log-log regression
- Power regression
- Exponential regression
- Percent logarithm
- Cubic equation regression
- 4 parameters regression**
- Quadratic regression

Double standard

Qualitative way

Negative < X

Positive >= X

Setting

New

Save

Delete

3.10 選擇樣本重複數

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Standard	Concentration

Double standard

Qualitative way

Negative < X

Positive >= X

Setting

New

Save

Delete

Item setting

QC setting

Test 450

Program STD 630

Full name

Unit ng/ml

wavelength 630 405

Reference 0 -- 1

Reagent Setting

Measurement End point meth

Calculation 4 parameters re

Samples 3

Intervals time

Detection

Factor

Standard setting

3.1.1 設定Standard

Standard number

0

Standard

Concentration

 Double standard**3.1.1.1 如Standard為兩重複請勾選**

Qualitative way None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

-
-
-
-
-
-
-
-
-

3.11.2 選擇Standard 1

Standard	Concentration
Standard1	

Qualitative way

Negative < X

Positive > = X

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

Standard	Concentration
Standard1	31.25

Qualitative way

3.11.3 輸入Standard 1濃度 (最低濃度)

e<

Positive>=

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

- 0
- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8

3.11.4 選擇Standard 2

Standard	Concentration
Standard1	31.25
Standard2	

Qualitative way

Negative < X

Positive >= X

Setting

New

Save

Delete

Item setting

QC setting

Test 450

Program

Full name

Unit

wavelength

Reference --

Reagent

Measurement

Calculation

Samples

Intervals time

Detection

Factor

Standard setting

Standard number

Double standard

Standard	Concentration
Standard1	31.25
Standard2	62.5

3.11.5 輸入Standard 2濃度

Qualitative way

Negative < X

Positive >= X

Item setting

QC setting

Test 450

Program STD 630

Full name

Unit ng/ml

wavelength 630 405

Reference 0 -- 1

Reagent

Setting

Measurement End point meth

Calculation 4 parameters re

Samples 3

Intervals time

Detection

Factor

Standard setting

Standard number 8

 Double standard

Standard	Concentration
Standard4	250
Standard5	500
Standard6	1000
Standard7	2000
Standard8	4000

Qualitative way None

Negative <

X

Positive >=

X

Setting

3.11.6 以此類推直到Standard 8設定完畢
(Standard 1為最低濃度, Standard 8為最高濃度)

Item setting

QC setting

Test 450

Program STD 630

Full name

Unit ng/ml

wavelength 630 405

Reference 0 -- 1

Reagent

Setting

Measurement End point meth

Calculation 4 parameters re

Samples 3

Intervals time

Detection

Factor

Standard setting

Standard number 8

 Double standard

Standard	Concentration
Standard4	250
Standard5	500
Standard6	1000
Standard7	2000
Standard8	4000

3.12 選擇定性方式（非定性分析則選None）

Qualitative way

None

None

Negative <

Positive threshold

Positive > =

X

Setting

New

Save

Delete

Item setting

QC setting

Test 450

Program Full name Unit wavelength Reference -- Reagent

Setting

Measurement Calculation Samples Intervals time Detection Factor

Standard setting

Standard number Double standard

Standard	Concentration
Standard4	250
Standard5	500
Standard6	1000
Standard7	2000
Standard8	4000

Qualitative way Negative < X Positive > = X

Setting

3.13 Program設定完成，按Save儲存。

New

Save

Delete



Item setting

QC setting

Test 450

STD 630

Program STD 630

Full name

Unit ng/ml

wavelength 630 405

Reference 0 -- 1

Reagent

Setting

3.14 儲存成功的Program會顯示在左側列表

Measurement End point meth

Calculation 4 parameters re

Samples 3

Intervals time

Detection

Factor

Standard setting

Standard number 8

 Double standard

Standard	Concentration
Standard2	62.5
Standard3	125
Standard4	250
Standard5	500
Standard6	1000

Qualitative way None

Negative <

X

Positive >=

X

Setting

New

Save

Delete

3.15 切換到Test介面

Item parameters

Test 450 **STD 630**



3.16 點選要套用的Program

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

Test parameters

Test template Setting Testing process Setting

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Item parameters

Test 450

STD 630

3.17 點選STD

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

3.18 選擇STD排列方向

Plate Direction

Row

Column

Start sampleNo

0001

Well Position

A

1

--

H

12

Setting

Plate current well position item information

3.19 依序點選空格以排列Standard位置

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1										
B	STD-2	STD-2										
C	STD-3	STD-3										
D	STD-4	STD-4										
E	STD-5	STD-5										
F	STD-6	STD-6										
G	STD-7	STD-7										
H	STD-8	STD-8										

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450

STD 630

3.20 點選Sample

Test parameters

Test template

Setting

Testing process

default

Setting

Sample NC PC PC2 CR QC1 QC2 QC3 QC4 QC5

STD Blank Clear

3.21 點選空格以排列Sample位置

Plate Direction Row Column

Start sampleNo 0001

Well Position A 1 -- H 12 Setting

Plate current well position item information

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001							
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							
H	STD-8	STD-8	0006	0007	0008							

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450

STD 630

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank

Clear

3.23 點選空格以排列Blank位置

3.22 點選Blank

Plate Direction

Row

Column

Start sampleNo

0001

Well Position

A

1

H

12

Setting

Plate current well position item information

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001	BK-1	BK-2					
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							
H	STD-8	STD-8	0006	0007	0008							

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

3.24 按Plate_Reset 兩次

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

Setting

Testing process

default

Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001	BK-1	BK-2					
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							
H	STD-8	STD-8	0006	0007	0008							

Item parameters

Test 450 STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

3.25 按Plate_Pop, 托架彈出, 放入盤子 (注意: AI朝右邊內側)。

Test parameters

Test template Testing process

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001	BK-1	BK-2					
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							
H	STD-8	STD-8	0006	0007	0008							

<input type="button" value="Plate_Pop"/>	<input type="button" value="Plate_Reset"/>	<input type="button" value="Test"/>	<input type="button" value="Stop testing"/>	<input type="button" value="Clear All"/>	<input type="button" value="Result"/>	<input type="button" value="ABS"/>	<input type="button" value="Quantitative"/>	<input type="button" value="Qualitative"/>
<input type="button" value="Save"/>	<input type="button" value="Print"/>	<input type="button" value="Print Preview"/>	<input type="button" value="Export"/>	<input type="button" value="Send"/>				



Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

3.26 按Test開始測讀

Test parameters

Test template Setting Testing process default Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001	BK-1	BK-2					
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							
H	STD-8	STD-8	0006	0007	0008							

Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

Test parameters

Test template

Testing process

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1	STD-1	0001	0001	0001	BK-1	BK-2					
B	STD-2	STD-2	0002	0002	0002							
C	STD-3	STD-3	0003	0003	0003							
D	STD-4	STD-4	0004	0004	0004							
E	STD-5	STD-5	0005	0005	0005							
F	STD-6	STD-6	0006	0007	0008							
G	STD-7	STD-7	0006	0007	0008							

(測讀過程中若要中止可按Stop testing)

<input type="button" value="Plate_Pop"/>	<input type="button" value="Plate_Reset"/>	<input type="button" value="Test"/>	<input type="button" value="Stop testing"/>	<input type="button" value="Clear All"/>	<input type="button" value="Result"/>	<input type="button" value="ABS"/>	<input type="button" value="Quantitative"/>	<input type="button" value="Qualitative"/>
<input type="button" value="Save"/>	<input type="button" value="Print"/>	<input type="button" value="Print Preview"/>	<input type="button" value="Export"/>	<input type="button" value="Send"/>				

Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo 0001

Well Position A 1 -- H 12 Setting

Plate current well position item information

Empty text box for plate information

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

Test parameters

Test template

Setting

Testing process

default

Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250	0.000	0.000					
B	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
C	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
D	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
E	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
F	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
G	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
H	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							

3.27 測讀完成，顯示結果。
(同一格內上為吸光值下為濃度)

Item parameters

Test 450 STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information


Test parameters

Test template Testing process

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250	0.000	0.000					
B	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
C	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
D	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
E	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
F	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
G	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							
H	0.000	0.000	0.000 <31.250	0.000 <31.250	0.000 <31.250							

<input type="button" value="Plate_Pop"/>	<input type="button" value="Plate_Reset"/>	<input type="button" value="Test"/>	<input type="button" value="Stop testing"/>	<input type="button" value="Clear All"/>	<input type="button" value="Result"/>	<input type="button" value="ABS"/>	<input type="button" value="Quantitative"/>	<input type="button" value="Qualitative"/>
<input type="button" value="Save"/>	<input type="button" value="Print"/>	<input type="button" value="Print Preview"/>	<input type="button" value="Export"/>	<input type="button" value="Send"/>				



3.28 按Result同時顯示吸光值及濃度 4:04:06

Item parameters

Test 450

STD 630

Plate Direction Row

Column

Start sampleNo

Well Position

A

1

--

H

12

Setting

Plate current well position item information

Empty text box for plate current well position item information.

Test parameters

Test template

Setting

Testing process

default

Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A	STD-1 0.000	STD-1 0.000	0001 0.000	0001 0.000	0001 0.000	BK-1 0.000	BK-2 0.000					
B	STD-2 0.000	STD-2 0.000	0002 0.000	0002 0.000	0002 0.000							
C	STD-3 0.000	STD-3 0.000	0003 0.000	0003 0.000	0003 0.000							
D	STD-4 0.000	STD-4 0.000	0004 0.000	0004 0.000	0004 0.000							
E	STD-5 0.000	STD-5 0.000	0005 0.000	0005 0.000	0005 0.000							
F	STD-6 0.000	STD-6 0.000	0006 0.000	0007 0.000	0008 0.000							
G	STD-7 0.000	STD-7 0.000	0006 0.000	0007 0.000	0008 0.000							
H	STD-8 0.000	STD-8 0.000	0006 0.000	0007 0.000	0008 0.000							

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

3.29 按ABS 顯示吸光值



Item parameters

Test 450

STD 630

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

Test parameters

Test template Setting Testing process Setting

Sample	NC	PC	PC2	CR	QC1	QC2	QC3	QC4	QC5
STD	Blank	Clear							

	1	2	3	4	5	6	7	8	9	10	11	12
A			0001 <31.250	0001 <31.250	0001 <31.250							
B			0002 <31.250	0002 <31.250	0002 <31.250							
C			0003 <31.250	0003 <31.250	0003 <31.250							
D			0004 <31.250	0004 <31.250	0004 <31.250							
E			0005 <31.250	0005 <31.250	0005 <31.250							
F			0006 <31.250	0007 <31.250	0008 <31.250							
G			0006 <31.250	0007 <31.250	0008 <31.250							
H			0006 <31.250	0007 <31.250	0008 <31.250							

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

Save

Print

Print Preview

Export

Send

3.30 按Quantitative 顯示濃度



Item parameters

Test 450

STD 630

3.3 | 按Save存档 (一定要先Save才能輸出CSV档)

Plate Direction Row Column

Start sampleNo

Well Position --

Plate current well position item information

Empty text box for plate information

Plate_Pop

Plate_Reset

Test

Stop testing

Clear All

Result

ABS

Quantitative

Qualitative

1

Save

Print

Print Preview

Export

Send

Test parameters

Test template

Setting

Testing process

default

Setting

Sample

NC

PC

PC2

CR

QC1

QC2

QC3

QC4

QC5

STD

Blank


Clear

	1	2	3	4	5	6	7	8	9	10	11	12
A			0001 <31.250	0001 <31.250	0001 <31.250							
B			0002	0002	0002 <31.250							
C					0003 <31.250							
D					0004 <31.250							
E					0005 <31.250							
F			0006 <31.250	0007 <31.250	0008 <31.250							
G			0006 <31.250	0007 <31.250	0008 <31.250							
H			0006 <31.250	0007 <31.250	0008 <31.250							

Set Board Number

Plate 2

3 Save



3.32 按Export輸出成CSV檔


TEST1 3

Method End point method

Calculate Percent logarithm

ABS	Concentration
0.100	1.000
0.200	2.000
0.300	3.000
0.400	4.000

3.33 切換到Report介面查看並輸出Standard Curve



Fit the equation: $\lg(y) = A*(x/xMax)^3 + B*(x/xMax)^2 + C*(x/xMax) + D$, correlation coefficient = 1.000

4 按Print列印成PDF档

Print

Send

注意！ 每一個Parameter (Program/Protocol) 只會存档第一次test的Standard Curve, 如欲Print每一次定量test的Standard Curve, 則每一次test前都必須重新建立新的Parameter (Program/Protocol)。


3.34 Standard Curve輸出PDF報告範例

Test 02 Standard curve line

Measure Method: End point method

Calc. Method: Linear regression

Fit equation: $y = A * x + B$, $R=0.9991$



ABS	Concentration
0.0047	0.0000
0.2458	0.0500
0.5355	0.1000
1.0303	0.2000

4. 查看歷史檢測數據

Sample record

Plate record

STD inquiry

QC inquiry

4.1 查看歷史檢測數據一般操作

Test date: 2022-06-06 2

Measure method:

Item No.:

Tester:

Sender:

Auditor:

Temperature: 0.0 °C

Humidity: 0.0 %

Sample No	Item name	Item full name	ABS	Qualitative	Quantitative	U	No.	Plate No	Test Time
							1	2206060001	2022-06-06 11:35

- ① 切換到Report介面
- ② 切換到Plate record分頁
- ③ 選擇欲篩選的檢測日期：今天（Today Sample）、所有日期（All Samples）、特定日期範圍（Search Result）

Export Save Delete Print Print Preview Send Curve

Sample list 3

Today Sample Today Sample All Samples Search Result

Query

Sample record | Plate record | **STD inquiry** | QC inquiry

Test date: 2022-04-29

Measure method: End point method

Item No.: 1

Tester: Admin

Sender: [Dropdown]

Auditor: [Dropdown]

Temperature: 0.0 °C

Humidity: 0.0 %

Sample No	Item name	Item full name	ABS	Qualitative	Quantitative	U
STD-1	STD 492		0.0032		0.0000	
STD-1	STD 492		0.0033		0.0000	
0001	STD 492		0.0999		0.0962	
0001	STD 492		0.1011		0.0973	
STD-2	STD 492		0.2442		0.2500	
STD-2	STD 492		0.2440		0.2500	
STD-3	STD 492		0.5339		0.5000	
STD-3	STD 492		0.5334		0.5000	
STD-4	STD 492		1.0309		1.0000	
STD-4	STD 492		1.0265		1.0000	

⑥ 該盤內所有檢測數據會列在這裡

④ 篩出所有檢測盤紀錄

No.	Plate No	Test Time
1	2206060001	2022-06-06 11:3
2	2204290002	2022-04-29 10:2
3	2204290001	2022-04-29 10:1
4	2204280002	2022-04-28 18:0
5	2204280001	2022-04-28 18:0
6	2111110006	2021-11-11 18:4
7	2111110005	2021-11-11 18:3
8	2111110004	2021-11-11 18:2
9	2111110003	2021-11-11 18:1
10	2111110002	2021-11-11 18:1
11	2111110001	2021-11-11 10:4
12	2111100001	2021-11-10 18:4
13	2009280003	2020-09-28 15:3
14	2009280002	2020-09-28 14:5
15	2009280001	2020-09-28 14:3
16	202009160002	2020-09-16 16:0



⑤ 點選欲檢視的盤子

Export | Save | Delete | Print | Print Preview | Send | Curve | Sample list | **All Samples** | Query

⑦ 可針對該盤的數據進行其他操作，如：Export CSV、Print PDF等。

③ 如：選擇所有日期

Sample record

Plate record

STD inquiry

QC inquiry

4.2 查看特定日期範圍的歷史檢測數據

Test date: 2022-04-29

Measure method: End point method

Item No.: 1

Tester: Admin

Sender: [Dropdown]

Auditor: [Dropdown]

Temperature: 0.0 °C

Humidity: 0.0 %

Sample No	Item name	Item full name	ABS	Qualitative	Quantitative	Unit	No.	Plate No	Test Time
STD-1	STD 492		0.0032		0.0000		1	2206060001	2022-06-06 11:3
STD-1	STD 492		0.0033		0.0000		2	2204290002	2022-04-29 10:2
0001	STD 492		0.0999		0.0962		3	2204290001	2022-04-29 10:1
0001	STD 492		0.1011		0.0973		4	2204280002	2022-04-28 18:0
STD-2	S						5	2204280001	2022-04-28 18:0
STD-2	S						6	2111110006	2021-11-11 18:4
STD-3	S						7	2111110005	2021-11-11 18:3
STD-3	S						8	2111110004	2021-11-11 18:2
STD-4	S						9	2111110003	2021-11-11 18:1
STD-4	S						10	2111110002	2021-11-11 18:1
							11	2111110001	2021-11-11 10:4
							12	2111100001	2021-11-10 18:4
							13	2009280003	2020-09-28 15:3
							14	2009280002	2020-09-28 14:5
							15	2009280001	2020-09-28 14:3
							16	202009160002	2020-09-16 16:0

Query

Plate No: [Text Box]

Measure method: [Dropdown]

Test date: 2022/ 4/ 1 - 2022/ 4/30

Tester: [Dropdown]

Sender: [Dropdown]

Auditor: [Dropdown]

Ok Cancel

- 選擇Search Result
- 按Query
- 設定篩選的起迄日期
- 如有需要可以選擇特定操作者
- 按OK

Export Save Delete Print Print Preview Send Curve

Sample list: Search Result [Dropdown] Query

Sample record Plate record **STD inquiry** QC inquiry

Test date: 2022-04-28

Measure method: End point method

Item No.: 1

Tester: Admin

Sender:

Auditor:

Temperature: 0.0 °C


Humidity: 0.0 %

Sample No	Item name	Item full name	ABS	Qualitative	Quantitative	Unit
STD-1	S492_0428		0.0036		0.0000	
STD-1	S492_0428		0.0038		0.0000	
0002	S492_0428		0.1052		0.1009	
0001	S492_0428		0.1058		0.1015	
STD-2	S492_0428		0.2446		0.2500	
STD-2	S492_0428		0.2444		0.2500	
STD-3	S492_0428		0.5341		0.5000	
STD-3	S492_0428		0.5333		0.5000	
STD-4	S492_0428		1.0315		1.0000	
STD-4	S492_0428		1.0264		1.0000	

⑧ 該盤內所有檢測數據會列在這裡

⑥ 篩出該日期區間內所有檢測盤紀錄

No.	Plate No	Test Time
1	2204290002	2022-04-29 10:26:
2	2204290001	2022-04-29 10:15:
3	2204280002	2022-04-28 18:06:
4	2204280001	2022-04-28 18:04:



⑦ 點選欲檢視的盤子

Export Save Delete Print Print Preview Send Curve

Sample list Search Result Query

⑨ 可針對該盤的數據進行其他操作，如：Export CSV、Print PDF等。

5. 設定使用者帳號密碼

No.	UserName	User level
1	Admin	Administrator

UserName

User level

Enter password

Enter password again

2
New

3
Save

4
Delete

6. 查看操作紀錄

History

Workload

2

3

6.1 查看這台機器所有操作紀錄

Date 2022-06-06 -- 2022-06-06 No. All

Item All Result All

Statistics

Print

Print Preview

Export

4

Date	No.	Item	Result
2022-06-06 11:35:25	0001	Test 02	
2022-06-06 11:35:25	0002	Test 02	

5

Statistical results

Total 2

Suspicious 0

Positive 0

Abnormal Ratio 0%

6.2 查看特定操作者的操作紀錄

1

2

Date 2022-01-01 -- 2022-06-06

3

User name Admin

4

Statistics Print Print Preview Export

User name	Item	Samples
Admin	STD 492	3
Admin	S492_0428	2
Admin	A450	96
Admin	Test 02	2

5

Statistical results

Plates 5

Items 4

Samples 100