


Lenti-p53 siRNA Virus, High Titer

Cat.No.	Unit
LV619	2 x 100 µl

Cat. No.	LV619
Name	Lenti-p53 siRNA Virus, High Titer
Unit	2 x 100 µl
Description	High Titer (10 ⁹ IU/ml) Recombinant Lentivirus expressing p53 siRNA.
Application	Cell immortalization.
Expression System Type	Lentivirus
Vector Map	Lenti-p53 siRNA
Vector Map	



Caution This product is distributed for laboratory research only.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).